

MAINTENANCE AND CLEANING SAFE WRAP RANGE

It is compatible with any cleaning protocols and resistant to most chemical agents, alcohol, diluted acids, oils.

The product can be cleaned/disinfected by all conventional cleaning methods, using non-abrasive accessories, cleaning products, detergents or products currently used in healthcare environments.

No decrease of the activity after 365 cleaning operations with water, ethanol, bleach or «ANIOSURF».

For best results!



WARM
WATER



SOFT
DETERGENT



COTTON
CLOTH

Regular cleaning (weekly) is recommended to keep your product in a nice condition.

Use warm water with a soft detergent and cotton cloth.