

Ready-to-use disinfectant and cleaning product. For surfaces, medical fittings and equipment.

(in accordance with Directive 93/42/EEC: MDD)

Effective against Coronavirus, Mers, Sars & Influenza



Delivery forms / Packagung units

Unit sizes Flask 500 ml Packaging units
12 Flasks / Package

Art. no BACDL500

Product

Ready-to-use disinfectant product providing long term efficacy. Does not contain Aldehyde or Phenol.

Uses

For the disinfection and cleaning of "medical" areas in accordance with Directive 93/42/EEC (Medical Devices) and all types of surfaces in hospitals, doctor's practices, rehabilitation centres and retirement homes. Especially useful in areas demanding effective and long-lasting hygiene.

Particularly suitable for areas where unpleasant odours caused by micro-organisms form, such as toilets and sanitary facilities. BacobanDL may be used in critical and sensitive areas of the pharmaceutical and cosmetic industries.

Microbiological effectiveness

Bacoban®DL is effective against: bacteria, fungi, viruses (hepatitis B and C, HIV, influenza including H5N1and H1N1, rotaviruses and adenoviruses)

Composition

100 g solution contains:

benzalkonium chloride 0,26 g, sodium pyrithione 0,025g, polycondenstaes, perfume substances, purified water.



| Effectiveness | |
|---------------|--------|
| DIN EN 1040 | 5 min. |
| DIN EN 1275 | 5 min. |
| DIN EN 1276 | 5 min. |
| DIN EN 1650 | 5 min. |
| DIN EN 13697 | 5 min. |

Long lasting anti-viral and anti-microbial performance. Effective for up to 10 days. Complies with ASTM E 2180 standards.

Bacoban® DL is effective against: bacteria, fungi, viruses (hepatitis B and C, HIV, influenza including H5N1and H1N1, rotaviruses and adenoviruses)

Application

Ready-to-use solution. To ensure optimum disinfectant performance, ensure that the whole area is coated. After use, thoroughly rinse cleaning aids with water. Do not mix with other cleaners or additives. Suitable for use on most surfaces. If unsure about the resistance of a particular surface, we recommend testing the product on an inconspicuous area. In some instances a visible coating may develop after drying, eg. when used on transparent surfaces or when over applied. If required this residue can be easily removed with standard commercial acidic or alkaline cleaners.

Always store in closed container. Keep cool and dry.

Container may not be recycled until completely empty.

Do not use after the expiry date. Dispose of correctly, in accordance with local and official regulations.

Safety advice:

Throughout the world legislation stipulates that appropriate protective measures and clothing be utilized when handling hazardous substances, therefore appropriate consideration should be made when applying Bacoban. In the event of contact with skin or eyes, rinse immediately and thoroughly with water and contact a doctor. Do not swallow. Keep out of the reach of children.

Important!

Bacoban offers proven long-lasting anti-microbial and anti-viral performance; however, it is important to stress that continued and appropriate disinfection procedures should be maintained in all instances.

Use disinfection products carefully. Ensure that the application instructions are adhered to.

UBA No: 5704 0031 Biocide reg. No: N-34071

(E₀₄₈₁